

Standardised Setups (minimum equipment = **Bold** marked)





Standard	OEM	Scratch figure acc. to standard	Order No.	Product Description
GME 60280 Method A	General Motors / Opel	<p>GME 60280 (Verfahren A)</p>	<b>03240131</b>	<b>SCRATCH HARDNESS TESTER 430 P-I-Smart</b>
				alternatively
			03240231	SCRATCH HARDNESS TESTER 430 P-II-Smart
				alternatively
			03240331	SCRATCH HARDNESS TESTER 430 P-III-Smart
				alternatively
			03240431	SCRATCH HARDNESS TESTER 430 P-IV-Smart
	<b>05390332</b>	<b>Test tip acc. to ISO 1518-1 (Ø 1 mm)</b>		
	<b>05670132</b>	<b>Set of weights for low load range (6 pcs.)</b>		
GME 60280 Method B	General Motors / Opel	<p>GME 60280 (Verfahren B)</p>	<b>03240131</b>	<b>SCRATCH HARDNESS TESTER 430 P-I-Smart</b>
				alternatively
			03240231	SCRATCH HARDNESS TESTER 430 P-II-Smart
				alternatively
			03240331	SCRATCH HARDNESS TESTER 430 P-III-Smart
				alternatively
			03240431	SCRATCH HARDNESS TESTER 430 P-IV-Smart
	<b>05390532</b>	<b>Testing device for writing effect acc. to Opel</b>		
	<b>05670132</b>	<b>Set of weights for low load range (6 pcs.)</b>		
GMW 14688 Method A	General Motors / Opel	<p>GMW 14688 (Verfahren A)</p>	<b>03240131</b>	<b>SCRATCH HARDNESS TESTER 430 P-I-Smart</b>
				alternatively
			03240231	SCRATCH HARDNESS TESTER 430 P-II-Smart
				alternatively
			03240331	SCRATCH HARDNESS TESTER 430 P-III-Smart
				alternatively
			03240431	SCRATCH HARDNESS TESTER 430 P-IV-Smart
	<b>05390332</b>	<b>Test tip acc. to ISO 1518-1 (Ø 1 mm)</b>		
	<b>05670132</b>	<b>Set of weights for low load range (6 pcs.)</b>		
GMW 14688 Method B	General Motors / Opel	<p>GMW 14688 (Verfahren B)</p>	<b>03240131</b>	<b>SCRATCH HARDNESS TESTER 430 P-I-Smart</b>
				alternatively
			03240231	SCRATCH HARDNESS TESTER 430 P-II-Smart
				alternatively
			03240331	SCRATCH HARDNESS TESTER 430 P-III-Smart
				alternatively
			03240431	SCRATCH HARDNESS TESTER 430 P-IV-Smart
	<b>05390532</b>	<b>Testing device for writing effect acc. to Opel</b>		
	<b>05670132</b>	<b>Set of weights for low load range (6 pcs.)</b>		

Standardised Setups (minimum equipment = **Bold** marked)

Standard	OEM	Scratch figure acc. to standard	Order No.	Product Description
PV 3952	Volkswagen / Audi		<b>03240131</b>	<b>SCRATCH HARDNESS TESTER 430 P-I-Smart</b>
				alternatively
			03240231	SCRATCH HARDNESS TESTER 430 P-II-Smart
				alternatively
			03240331	SCRATCH HARDNESS TESTER 430 P-III-Smart
				alternatively
			03240431	SCRATCH HARDNESS TESTER 430 P-IV-Smart
	<b>05390332</b>	<b>Test tip acc. to ISO 1518-1 (Ø 1 mm)</b>		
	<b>05670132</b>	<b>Set of weights for low load range (6 pcs.)</b>		
D45 1047	Renault		<b>03240131</b>	<b>SCRATCH HARDNESS TESTER 430 P-I-Smart</b>
				alternatively
			03240231	SCRATCH HARDNESS TESTER 430 P-II-Smart
				alternatively
			03240331	SCRATCH HARDNESS TESTER 430 P-III-Smart
				alternatively
			03240431	SCRATCH HARDNESS TESTER 430 P-IV-Smart
	<b>05390332</b>	<b>Test tip acc. to ISO 1518-1 (Ø 1 mm)</b>		
	<b>05670132</b>	<b>Set of weights for low load range (6 pcs.)</b>		
MS210-05K /4.9.1	Hyundai / Kia		<b>03240131</b>	<b>SCRATCH HARDNESS TESTER 430 P-I-Smart</b>
				alternatively
			03240231	SCRATCH HARDNESS TESTER 430 P-II-Smart
				alternatively
			03240331	SCRATCH HARDNESS TESTER 430 P-III-Smart
				alternatively
			03240431	SCRATCH HARDNESS TESTER 430 P-IV-Smart
	<b>05390332</b>	<b>Test tip acc. to ISO 1518-1 (Ø 1 mm)</b>		
	<b>05670132</b>	<b>Set of weights for low load range (6 pcs.)</b>		
MS210-06K	Hyundai / Kia		<b>03240131</b>	<b>SCRATCH HARDNESS TESTER 430 P-I-Smart</b>
				alternatively
			03240231	SCRATCH HARDNESS TESTER 430 P-II-Smart
				alternatively
			03240331	SCRATCH HARDNESS TESTER 430 P-III-Smart
				alternatively
			03240431	SCRATCH HARDNESS TESTER 430 P-IV-Smart
	<b>05390332</b>	<b>Test tip acc. to ISO 1518-1 (Ø 1 mm)</b>		
	<b>05670132</b>	<b>Set of weights for low load range (6 pcs.)</b>		

Standardised Setups (minimum equipment = **Bold** marked)

Standard	OEM	Scratch figure acc. to standard	Order No.	Product Description
STD 4377	Volvo	<p>STD 4377</p>	<b>03240131</b>	<b>SCRATCH HARDNESS TESTER 430 P-I-Smart</b>
				alternatively
			03240231	SCRATCH HARDNESS TESTER 430 P-II-Smart
				alternatively
			03240331	SCRATCH HARDNESS TESTER 430 P-III-Smart
				alternatively
			03240431	SCRATCH HARDNESS TESTER 430 P-IV-Smart
		<b>05390332</b>	<b>Test tip acc. to ISO 1518-1 (Ø 1 mm)</b>	
		<b>05670132</b>	<b>Set of weights for low load range (6 pcs.)</b>	
PV 3974	Volkswagen / Audi	<p>PV 3974</p>	<b>03240131</b>	<b>SCRATCH HARDNESS TESTER 430 P-I-Smart</b>
				alternatively
			03240231	SCRATCH HARDNESS TESTER 430 P-II-Smart
				alternatively
			03240331	SCRATCH HARDNESS TESTER 430 P-III-Smart
				alternatively
			03240431	SCRATCH HARDNESS TESTER 430 P-IV-Smart
		<b>05391132</b>	<b>Testing device for writing effect acc. to VW/Audi</b>	
		<b>05670132</b>	<b>Set of weights for low load range (6 pcs.)</b>	
PV 3962	Volkswagen / Audi	<p>PV 3962</p>	<b>03240131</b>	<b>SCRATCH HARDNESS TESTER 430 P-I-Smart</b>
				alternatively
			03240231	SCRATCH HARDNESS TESTER 430 P-II-Smart
				alternatively
			03240331	SCRATCH HARDNESS TESTER 430 P-III-Smart
				alternatively
			03240431	SCRATCH HARDNESS TESTER 430 P-IV-Smart
		<b>05391432</b>	<b>Test tip acc. to PV 3962 (VW/Audi)</b>	
		<b>05670132</b>	<b>Set of weights for low load range (6 pcs.)</b>	
DBL 9202	Daimler	<p>DBL 9202</p>	<b>03240131</b>	<b>SCRATCH HARDNESS TESTER 430 P-I-Smart</b>
				alternatively
			03240231	SCRATCH HARDNESS TESTER 430 P-II-Smart
				alternatively
			03240331	SCRATCH HARDNESS TESTER 430 P-III-Smart
				alternatively
			03240431	SCRATCH HARDNESS TESTER 430 P-IV-Smart
		<b>05390932</b>	<b>Testing device for test acc. to DBL 9202</b>	
		<b>05670132</b>	<b>Set of weights for low load range (6 pcs.)</b>	

Standardised Setups (minimum equipment = <b>Bold</b> marked)				
Standard	OEM	Scratch figure acc. to standard	Order No.	Product Description
GS 97034-8	BMW	<p>GS 97034-8</p> <p>F [N]</p> <p>1 _____</p> <p>3 _____</p> <p>5 _____</p> <p>8 _____</p> <p>10 _____</p> <p>12 _____</p> <p>15 _____</p> <p>20 _____</p> 	<b>03240131</b>	<b>SCRATCH HARDNESS TESTER 430 P-I-Smart</b>
				alternatively
			03240231	SCRATCH HARDNESS TESTER 430 P-II-Smart
				alternatively
			03240331	SCRATCH HARDNESS TESTER 430 P-III-Smart
				alternatively
			03240431	SCRATCH HARDNESS TESTER 430 P-IV-Smart
			<b>05390432</b>	<b>Test tip acc. to BMW (Ø 3.0 mm)</b>
			<b>05670132</b>	<b>Set of weights for low load range (6 pcs.)</b>
GS 97034-9	BMW	<p>GS 97034-9</p> <p>F [N]</p> <p>1 _____</p> <p>3 _____</p> <p>5 _____</p> <p>8 _____</p> <p>10 _____</p> <p>12 _____</p> <p>15 _____</p> <p>20 _____</p> 	<b>03240131</b>	<b>SCRATCH HARDNESS TESTER 430 P-I-Smart</b>
				alternatively
			03240231	SCRATCH HARDNESS TESTER 430 P-II-Smart
				alternatively
			03240331	SCRATCH HARDNESS TESTER 430 P-III-Smart
				alternatively
			03240431	SCRATCH HARDNESS TESTER 430 P-IV-Smart
			<b>05390132</b>	<b>Test tip acc. to van Laar (Ø 0.5 mm)</b>
			<b>05670132</b>	<b>Set of weights for low load range (6 pcs.)</b>
TPJLR.52.008	Jaguar/ Landrover	<p>TPJLR.52.008</p> 	<b>03240131</b>	<b>SCRATCH HARDNESS TESTER 430 P-I-Smart</b>
				alternatively
			03240231	SCRATCH HARDNESS TESTER 430 P-II-Smart
				alternatively
			03240331	SCRATCH HARDNESS TESTER 430 P-III-Smart
				alternatively
			03240431	SCRATCH HARDNESS TESTER 430 P-IV-Smart
<b>05390132</b>	<b>Test tip acc. to van Laar (Ø 0.5 mm)</b>			
<b>05670132</b>	<b>Set of weights for low load range (6 pcs.)</b>			
TPJLR.52.010	Jaguar/ Landrover	<p>TPJLR.52.010</p> 	<b>03240131</b>	<b>SCRATCH HARDNESS TESTER 430 P-I-Smart</b>
				alternatively
			03240231	SCRATCH HARDNESS TESTER 430 P-II-Smart
				alternatively
			03240331	SCRATCH HARDNESS TESTER 430 P-III-Smart
				alternatively
			03240431	SCRATCH HARDNESS TESTER 430 P-IV-Smart
<b>05391232</b>	<b>Testing device for writing effect acc. to Jaguar/Landrover</b>			
<b>05670132</b>	<b>Set of weights for low load range (6 pcs.)</b>			